

ZARROSEN®

Desensitizing Cavity Varnish and Dentinal Tubuli Seal

Zarosen's unique, strontium-chloride formula provides immediate, long-lasting relief from hypersensitivity discomfort. Zarosen is effective for weeks when applied to the enamel, dentin or cementum of the tooth with the cotton pledget or microbrush. Use Zarosen after prophylaxis, scaling, tooth bleaching, restoration preparation, crown and bridge cementation and pin and post seating. Zarosen's special solvent system retains solids in solution and does not require the addition of a thinner to maintain original viscosity. Do not use Zarosen with resin-based cement or filling material.




Fig 1. Sensitive tooth



Fig 2. Swab Zarosen onto the tooth




Fig 3. Invisible barrier is now created providing long lasting protection from tooth hypersensitivity



Description

Zarosen® Desensitizing Cavity Varnish and Dentinal Tubuli Seal

Size

14 g

Item No.

0800

Key Features

- One bottle contains approximately 100 applications
- Fast-drying within 30 seconds
- One-step procedure
- Penetrates, fills and seals tubuli
- Tooth need not be dried prior to application
- May be applied to cementum, dentin, or enamel tooth structure
- No burnishing required

Benefits

- Won't interfere with patient eating or drinking after procedure
- Effectively reduces sensitivity of tooth bleaching procedures
- Lasts for weeks
- Indicated for a variety of procedures
- No thinner required to maintain viscosity
- Immediate relief

Cetylite Inc.
9051 River Road
Pennsauken, NJ 08110
856.665.6111
marketing@cetylite.com

©2022 Cetylite Industries, Inc. All rights reserved. Information is summary in nature and subject to change. Zarosen and Cetylite are registered trademarks of Cetylite Industries, Inc.

12CETY2022