

Cetylcide-II®

CONCENTRATE - Hard Surface Disinfectant

Our economically-priced Cetylcide-II Concentrate is an EPA-registered, broad-spectrum hospital disinfectant. It is germicidal, fungicidal and virucidal. Cetylcide-II is effective against 140+ disease-causing organisms including; HIV, MRSA, VRE, HBV, and Influenza A/Brazil Virus (H1N1).* Easily diluted in tap water, Cetylcide-II concentrate is now available in FOUR convenient packages to meet the needs of facilities of every size, from single operator to surgicenter to full size hospital.

Description	Yield	Item No.
32 oz Concentrate	64 Qts/16 Gal.	0152
12 Pack 1/2 oz bottles	12 Qts/3 Gal.	0159
6 Pack 1/2 oz bottles	6 Qts./1.5 Gal.	0158
1/2 oz bottle with Quart Spray Bottle	1 Quart	0153



Key Features

- Kills 140+ pathogenic organisms including HIV, Influenza A, Hepatitis B and C, VRE, MRSA & Influenza A/Brazil Virus (H1N1)*
- 1/2 oz bottle makes 1 quart of ready-to-use solution
- 2 oz of Cetylcide-II Concentrate makes 1 gallon of use-solution
- 1 quart of Cetylcide-II Concentrate makes 16 gallons of use solution
- No-mess dispensing with easy tip-and-pour bottle
- Broad spectrum surface disinfectant/deodorizer
- Complies with OSHA Bloodborne Pathogen Rule
- Pleasant lemon scent
- No alcohol or phenol

Benefits

- Economical, up to 80% less expensive than ready-to-use formulas
- Saves inventory space
- No waste - mix only what you need
- Additional Quart Spray Bottles may be purchased separately
- Dilutes completely in ordinary tap water

Contact:

Cetylite Inc.
9051 River Road
Pennsauken, NJ 08110
856.665.6111
marketing@cetylite.com

SARS-CoV-2 (COVID-19)

CDC recommends disinfection in healthcare settings using EPA-registered hospital-grade disinfectants shown on EPA's List N: Disinfectants for use against SARS-CoV-2. (Cetylcide-II Primary Product EPA Registration #61178-1)
<https://www.epa.gov/pesticide-registration/list-n-disinfectants-use-against-sars-cov-2>

*Except in California



	Cetyl-G® High Level Disinfectant / Sterilant & Diluent (item #0122)	Cetyl-II® Hard Surface Disinfectant (item #s 0152, 0153, 0158, 0159)	Cetylite Power Cleaner® Ultrasonic Cleaning Solution (item #0130)	Cetyl-Zyme Pro-Am® Dual Enzymatic Detergent Concentrate (item #0170)	Cetyl-Zyme Pro-Am® Dual Enzymatic Foam Spray (item #0176)	Hyde-Out® Aldehyde Neutralizer (item #0128)
Holding Solution			●	●		
Ultrasonic Solution			●	●		
Pre-Cleaning Solution			●	●	●	
Cleaning Solution		●	●			
Hard Surface Disinfectant		●				
Evacuation Line Cleaner				●		
High Level Disinfectant	●					
Sterilant	●					
Glutaraldehyde/OPA Neutralizer						●

Mixing Instructions

- 1 1/2 oz. of Cetyl-II Concentrate

OR
- 2 Fill Quart Spray Bottle with 1 quart of water.
- 3 Add 1/2 oz. (15 mL) of Cetyl-II Concentrate into Quart Spray Bottle.

OR
- 4 Secure trigger spray to bottle.

Refer to the product insert for complete directions for use.

©2020 Cetylite Industries, Inc. All rights reserved. Information is summary in nature and subject to change. Cetyl-Zyme Pro-Am, Cetyl-II, Cetyl-G, Hyde-Out, C-Tub, and Cetylite are registered trademarks of Cetylite Industries, Inc. All other copyrights are the property of their respective owners.

